



Customer Account Preferences

Default preferences are highlighted in BOLD

NOTE: WRITTEN INSTRUCTIONS on Rx will always override preferences

Dr. Name: _____ Cell # to text: _____

Address: _____

Email: _____ Website: _____

Phone #: _____ Office Hours: _____

For tech questions, who should we contact? Doctor Assistant Other: _____

How would you like us to contact you? email text phone

You can change your preferences at any time, by calling **(877) 452-4767**

of coats of die spacer? Circle one: 0 1 **2** 3 4




* Dr. to specify additional coats if desired.

Foil/Paint opposing? Yes **No**

What type margin do you typically prep?

Margin Prep (Select one): Butt Joint Champher
 Bevel Feather Edge

What is your preferred metal margin design?

Metal Design: Metal Lingual Collar 
 No Metal Collar 
 Metal Collar 360 

Our restorations are made to hold articulating paper when all other teeth are in occlusion. Is this acceptable?

Occlusion (Select one): In Occlusion **Slightly Out of Occlusion***
 Out of Occlusion Visibly Out of Occlusion

*Slightly out = One thickness of Madam Butterfly silk pulls with no resistance.

Pontic Design (Select one):



Hygienic



Bullet



Modified
Ridge Lap



Ridge Lap



Ovate

Our restorations are made to have a snap floss proximal contact. Is this acceptable?

Contacts (Select one): Broad and Light **Broad and Firm** Broad and Tight

Do you prefer occlusal pit and fissure staining on your PFM cases?

Pit and Fissure Staining (Select one):

Light = only stain the pits

Medium = stain the pits and a few secondary grooves

Heavy = stain the pits and secondary grooves

***None = only glaze the restoration**

Light

Medium

Heavy

None*

Do you currently mill crowns in your office? Yes No

If yes, what system do you use? _____

Do you have an IOS scanner? Yes No

If yes, what system do you have? _____

General Comments: _____

Please email this form to info@kairosdental.com or fax it to (714) 754-5452